

THYN1-thymocyte nuclear protein 1 Gene

Reference: PTXBC006978

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC144; MDS012; MY105; THY28; THY28KD; thymocyte nuclear protein 1; thymocyte protein thy28		

For research use only.

GENE DESCRIPTION	thymocyte nuclear protein 1
GENE ID	29087
ORIGIN SPECIES	Human
ACCESSIONS	BC006978
VECTOR	pDONR223
NCBI SYMBOL	THYN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgagaccccggaagaggctggctgggacttctgttcagacaagggaactatcaggaaaacgcacaaaactgagaactcaggtgaggcattagctaaagtggaggactccaaccctcagaagactcagccactaaaaactgtttgaagaatctaagcagccactggctgatgaa
gtcagagccagagagccgcttagagaaggtgtatgtgaaagttcagcattgaggtatctcaaagcacagccaaacagacaacatgctgggatgtgttcgtaactaccaggctcggaacttccttagagccatgaagctgggagaagaagccttcttaccatagcaactgcaaga
gccaggcatcgaggactcatgaagatcgtgaaagaggcttaccagaccacacagtttgagaaaaacaatccccattatgacccatctagcaaagaggacaaccctaagtgtccatgaagagtttgatttttttga