

CRIP1-cysteine-rich PDZ-binding protein Gene

Reference: PTXBC006980

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	postsynaptic protein CRIP1; HSPC139; SSMD1; cysteine-rich PDZ-binding protein; cysteine-rich interactor of PDZ three; cysteine-rich interactor of PDZ3; CXXC repeat containing interactor of PDZ3 domain		

For research use only.

GENE DESCRIPTION	cysteine-rich PDZ-binding protein
GENE ID	9419
ORIGIN SPECIES	Human
ACCESSIONS	BC006980
VECTOR	pDONR223
NCBI SYMBOL	CRIP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgtgcgaaaaatgtgaaagaaacttggtactgttatcactccagatcatggaagatggtgctaggaataccacagaaagtgtggaagaaagctgaatgaaataaagcttgactcaaaaaagcaagattgatccataggaagaataagttccactgtagaattgttaaagttctgtgcaccaaccagggttctcattactgccagggtgtgctcacaataaaggcatctgtgcgatgtgtggaataaaggttttgataccaaaactacaagcaaacatctgtctag