

TIMM8A-translocase of inner mitochondrial membrane 8 homolog A (yeast) Gene

Reference: PTXBC006994

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DDP; DDP1; DFN1; MTS; TIM8; mitochondrial import inner membrane translocase subunit Tim8 A; X-linked deafness dystonia protein; deafness dystonia protein 1; deafness/dystonia peptide; translocase of inner mitochondrial membrane 8 homolog A (yeast)		

For research use only.

GENE DESCRIPTION	translocase of inner mitochondrial membrane 8 homolog A (yeast)
GENE ID	1678
ORIGIN SPECIES	Human
ACCESSIONS	BC006994
VECTOR	pDONR223
NCBI SYMBOL	TIMM8A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattctctctcttccctccgcccgggttgggtgcagtgacccgcagttgcagcatctcatcgaggtagagactcaaaagcagcgttcacagcagctggtgcaccagatgactgaactttgttgggagaagtgcattggacaagcctggccaaagttggacagtcgggctgagg
cctgtttgtgaactgcgttgagcgttcattgatacaagccagttcatcttgaatcgactggaacagaccagaaatccaagccagttttctcagaagcctttctgactga