

VAMP4-vesicle-associated membrane protein 4 Gene

Reference: PTXBC007019

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VAMP-4; VAMP24; vesicle-associated membrane protein 4; vesicle associated membrane protein 4		

For research use only.

GENE DESCRIPTION	vesicle-associated membrane protein 4
GENE ID	8674
ORIGIN SPECIES	Human
ACCESSIONS	BC007019
VECTOR	pDONR223
NCBI SYMBOL	VAMP4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctcccaagttaagcgccacctcaatgatgatgtcacaggtctgtgaaaagtgaaggagaaatcttttgaagatgattcagatgaagaagaggacttttttggaccatctggaccaagatttggacctagaaatgataaaattaagcatgttcagaatcaagtggatgaagtt
attgatgtcatgcaagaaatattacaaggttaattgagagaggggagagactagatgaactacaggacaaatcagaagccttatcgataatgcaacagcttttagcaacagatccaacaacttcgaagcgaatgtgtgtggcgtggatgcaaaataaaagccatcatggctttgtgtgct
gctatcctttgtctagtattatcattcttatagtcataataccgtacttga