

SAA1-serum amyloid A1 Gene

Reference: PTXBC007022

| | | | |
|----------------|--|---------------------|----------------|
| Product type: | DNA | Storage Condition: | -20°C |
| Size/Quantity: | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms | PIG4; SAA; SAA2; TP5314; serum amyloid A-1 protein; serum amyloid A protein; tumor protein p53 inducible protein 4; serum amyloid A1 | | |

For research use only.

| | |
|-------------------------|---|
| GENE DESCRIPTION | serum amyloid A1 |
| GENE ID | 6288 |
| ORIGIN SPECIES | Human |
| ACCESSIONS | BC007022 |
| VECTOR | pDONR223 |
| NCBI SYMBOL | SAA1 |
| LENGTH | ProteoGenix bp |
| SEQUENCE PRIMERS | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

ORF (Only for reference, product provided without STOP codon)

atgaagcttctcacggcctggtttctgtccttggctcgtgggtgtcagcagccgaagctctttctccttggcagggctttgatgggctcgggacatgtggagcctactctgacatgagagaagccaattacatcggtcagacaatacttccatgctcggggaactatgatgctgccaaggggacctgggggtctgtggctgcagaagcgtcagcgtatgccagagagaatccagagattcttggccatgtgtcggaggactcactgccgatcagctgccgatgaatgggcaggagtggcaagacccaatcacttccacctgctggcctgcctgagaa tactga