

ELA3A-elastase 3A, pancreatic Gene

Reference: PTXBC007028

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ELA3A; ELA3; chymotrypsin-like elastase family member 3A; elastase 1; elastase 3A, pancreatic; elastase IIIA; elastase-3A; protease E; chymotrypsin like elastase family member 3A		

For research use only.

GENE DESCRIPTION	elastase 3A, pancreatic
GENE ID	10136
ORIGIN SPECIES	Human
ACCESSIONS	BC007028
VECTOR	pDONR223
NCBI SYMBOL	ELA3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgccggtgctcagttccctctcttctgtggcgttgctcaggtatggccaccttctctcactctccagcgcgttgccatggtgaggatgcggtccctacagctggccctggcaggttccctgcagatgagaaaagtgaagcttaccacacgtgtggcgtagcct
catcgcctccgaltgggttgactgccgcccactgcctcaggggatctgacctaccaggtggtgttgggtgagtaaaccttgcctgaaaggaggcccccagcaggtgatcccatcaactctgaggagctgtttgcatccactctggaaccgctcgtgtgtggcctgtggcaatg
acatgccctcatcaagctctcagcagcggccagctgggagatgccgtccagctcgcctcactccctccgctggtgacatctctcccaacaagacacctgctacatcaccgctggggcgtctctataccaatgggccactcccagacaagctgcagcagcccgctgccgtg
gtggactataagcactgctccaggtggaaactggtgggttccaccgtgaagaaaacatggtgtgtgctgagggtacatccgctccggtgcaacggtgactctggaggaccctcaactgccccacagaggatggtggctggcaggtccacgggtgtgaccagctttgttctgcttgc
gctgcaacttcatctggaagcctacagtggtcactcaggtctccgcttcatcactggtgaggagaccatagcaagccactag