

MT2A-metallothionein 2A Gene

Reference: PTXBC007034

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MT2; metallothionein-2; MT-2; MT-II; metallothionein-II; metallothionein 2A		

For research use only.

GENE DESCRIPTION	metallothionein 2A
GENE ID	4502
ORIGIN SPECIES	Human
ACCESSIONS	BC007034
VECTOR	pDONR223
NCBI SYMBOL	MT2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatcccaactgctcctgcgcgcgggtgactcctgcacctgcgcgcggctcctgcaaatgcaagagtgcaaatgcacctcctgcaagaaaagctgctcctgctgcctgctgcctgctggcctgtgccagtgtgccagggtgcacatgcaaggggcgtcgacaagtgcagctgctgcgcctga