



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CD69-CD69 molecule Gene

Reference: PTXBC007037

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD69 molecule; activation inducer molecule (AIM/CD69); CD69 antigen (p60, early T-cell activation antigen); early activation antigen CD69; AIM; BL-AC/P26; CLEC2C; EA1; GP32/28; MLR-3; C-type lectin domain family 2, member C; early T-cell activation antigen		

For research use only.

GENE DESCRIPTION	CD69 molecule
GENE ID	969
ORIGIN SPECIES	Human
ACCESSIONS	BC007037
VECTOR	pDONR223
NCBI SYMBOL	CD69
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr