

SLA-Src-like-adaptor Gene

Reference: PTXBC007042

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SLA1; SLAP; src-like-adaptor; SLAP-1; hSLAP; src-like-adaptor protein 1; Src-like-adaptor		

For research use only.

GENE DESCRIPTION	Src-like-adaptor
GENE ID	6503
ORIGIN SPECIES	Human
ACCESSIONS	BC007042
VECTOR	pDONR223
NCBI SYMBOL	SLA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaacagcatgaaatccacccctgcgcctgcccagagggccctgcccacccggaggactggatagcgaactccttgcctgctaagtgaactccctgacatcagcccccgatattccgaggggagaaactgcgtgtgattctgatgaaggggctggtgaa
agctatttctcttagcactggctcagagaggttacatccctggaatatggtgcccagagtttaccatggctgctgtttgagggcctgggagagacaagccgagggagctgctgcagctgccagacacaaaggctggctctttatgatcagagagagtgagaccaagaagggtttact
cactgtcgggtgagacacaggcaggtaaagcattaccgcatcttccgtctgcccacaactggtaactacatttccccgaggtcacttccagtgcctggaggacctggtgaaccattatctgaggtggctgatggcctgtgctgtgtcaccacggcctgcctgacacaaagcacggctg
ccccagcagtgagggcctccagctcactgtcacctgtcgtcagaagactgtgactggaggagagtgctcagactgcaggaggaccccgagggaacagagaacccgcttgggtgtagacgagtccttttcagctatggccttcgagagagcattgcctttacgttccctgaccagtg
aggacaacacctctttgatcgaagaagaaaaagcatctccctgatgtatgggtgagcaagagaagagctcatttctcatcaccacctactttgaggactag