

RGS2-regulator of G-protein signaling 2, 24kDa Gene

Reference: PTXBC007049

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	G0S8; regulator of G-protein signaling 2; G0 to G1 switch regulatory 8, 24kD; G0/G1 switch regulatory protein 8; cell growth-inhibiting gene 31 protein; cell growth-inhibiting protein 31; regulator of G-protein signaling 2, 24kDa		

For research use only.

GENE DESCRIPTION	regulator of G-protein signaling 2, 24kDa
GENE ID	5997
ORIGIN SPECIES	Human
ACCESSIONS	BC007049
VECTOR	pDONR223
NCBI SYMBOL	RGS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaaatgtctatgttcttgctgttcaacacgactgcagaccatggacaagagcgcaggcagtgccacagagcagagagaagcagagaaagatgaacggaccctttaaaagattggaagaccgtttgagctactcttcaaaaattcctactcctgggaagcccaaacg
gcaaaaaagcaaacagcaagctttcatcaagccttctcctgaggaagcagcgtgtgtcagaagcatttgacgagctgctagccagcaaatatgtcttgcctcattcaggccttttttaagtcggaattctgtgaagaaatattgaattctgctggcctgtgaagactcaaaaaacc
aatcaccaccaaaagctgtcctcaaaagcaaggaataatatactgacttcatagaaaaaggaagctcctcaaaagagataaacaatagattttcaaaccaaaactctgattgccagaatatacaagaagctacaagtggtgctttacaactgccagaaaagggtatagcttgatggagaac
aactcttctcctgtttcttgagtcagaattctaccaggactgtgttaaaagcccaaatcaccacagagcctcatgctacatga