

VAMP3-vesicle-associated membrane protein 3 (cellubrevin) Gene

Reference: PTXBC007050

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CEB; vesicle-associated membrane protein 3; VAMP-3; cellubrevin; synaptobrevin-3; vesicle associated membrane protein 3		

For research use only.

GENE DESCRIPTION	vesicle-associated membrane protein 3 (cellubrevin)
GENE ID	9341
ORIGIN SPECIES	Human
ACCESSIONS	BC007050
VECTOR	pDONR223
NCBI SYMBOL	VAMP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctacaggtccaactgctgccactggcagtaatcgaagacttcagcagacacaaatcaagtagatgaggtggtggacataatgcgagttacgtggacaaggttctggaaaagagaccagaagctctctgagttagacaccgtgcagacgcactgcaggcaggcgcttctcaattg
aaacgagcgcagccaagttgaagaggaaatattggtggaagaattgcaagatgtggcaatcgggattactgttctgttatcttcatcatcatcatcatctgtgtgggtgtctcttcattga