

UBL5-ubiquitin-like 5 Gene

Reference: PTXBC007053

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HUB1; ubiquitin-like protein 5; beacon; testicular tissue protein Li 217; ubiquitin like 5		

For research use only.

GENE DESCRIPTION	ubiquitin-like 5
GENE ID	59286
ORIGIN SPECIES	Human
ACCESSIONS	BC007053
VECTOR	pDONR223
NCBI SYMBOL	UBL5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcgagggtgttgcacgaccgtctggggaagaaggtccgctgtaaatgcaacacggatgataccatcgggacctaagaagctgattgcagccaaactggtaccgttggacaagattgtcctgaagaagtggtacacgattttaaggaccacgtgtctctggggactatga
aatccacgatgggatgaacctggagctttattatcaatag