

LSM3-LSM3 homolog, U6 small nuclear RNA associated (*S. cerevisiae*) Gene

Reference: PTXBC007055

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LSM3 homolog, U6 small nuclear RNA and mRNA degradation associated; LSM3 homolog, U6 small nuclear RNA associated; LSM3 U6 small nuclear RNA and mRNA degradation associated; U6 snRNA-associated Sm-like protein LSm3; SMX4; USS2; YLR438C		

For research use only.

GENE DESCRIPTION	LSM3 homolog, U6 small nuclear RNA associated (<i>S. cerevisiae</i>)
GENE ID	27258
ORIGIN SPECIES	Human
ACCESSIONS	BC007055
VECTOR	pDONR223
NCBI SYMBOL	LSM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggacgacgtagaccagcaacaactaccaacactgtagaggagccctggatcttatcaggctcagcctagatgagcgaatttatgtgaaaatgagaaatgaccgagagcttcgaggcagattacatgcttatgatcaacatttaaatatgatcttgggagatgtggaagaaactgtg
actactatagaattgatgaagaacatatgaagagatatataatcaacgaaacggaatatccaatgctctttgtccggggagatggcgttgtctgtgttgcctccctcactgagagttggctga