

ANKRD26L1-ankyrin repeat domain 26-like 1 Gene

Reference: PTXBC007072

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	ankyrin repeat domain 26-like 1
GENE ID	84832
ORIGIN SPECIES	Human
ACCESSIONS	BC007072
VECTOR	pDONR223
NCBI SYMBOL	ANKRD26L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaacaagaactcgcagtttgagatttctgatagtcataaaagaagaacgtgtgcataaaaccatttgatgcaagatgaaattgccaggctcaggctggaatacacacaataaaaaatcaaatcctgaaaagaatactaaaagacattgaaattataaaagaagcatg
aagacctcaaaaggctctaaacagaatggggaataatcaaaaaacgatagccattatagtgacagcttactgctctgacagatgagaacacaatgctcgttctaaactggagaagaaaaacaagcaggcaagactgacgaatggaatcataccatcgtagactga