

USMG5-up-regulated during skeletal muscle growth 5 homolog (mouse) Gene

Reference: PTXBC007087

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DAPIT; HCVFTP2; bA792D24.4; up-regulated during skeletal muscle growth protein 5; Diabetes Associated Protein in Insulin-sensitive Tissues; HCV F-transactivated protein 2; diabetes-associated protein in insulin-sensitive tissues; upregulated during skelet		

For research use only.

GENE DESCRIPTION	up-regulated during skeletal muscle growth 5 homolog (mouse)
GENE ID	84833
ORIGIN SPECIES	Human
ACCESSIONS	BC007087
VECTOR	pDONR223
NCBI SYMBOL	USMG5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaggtccagaaagtgatgcgaataccagttcactggtattaaaaatattcaactcttactctcacaggtagaatgaactgtgtactggccacataggagcattgcattgattgtcttatattcaagtttaaggtccaaaaaacccagctgtgaaagcaacataa