

ZNF626-zinc finger protein 626 Gene

Reference: PTXBC007116

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 626; CTC-513N18.7		
For research use only.			

GENE DESCRIPTION	zinc finger protein 626
GENE ID	199777
ORIGIN SPECIES	Human
ACCESSIONS	BC007116
VECTOR	pDONR223
NCBI SYMBOL	ZNF626
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggacattgcaatttagagatgtggccatagaattctctctggaggagtgccattgcctggacactgcacagcggaaatttatataggaatgtgatgttagagaactacagtaacctggctctctgtgtattactgtttctaagccagacctgatcacctgtctggagcaaggagaaaacctt
tgaccatgaagagaaatgagatgatagccaaaccctcagtgagcttcttcaagttcacagtgagagccaaagtctcttcatgacatataa