



Reference: PTXBC007213

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	palmitoyltransferase ZDHHC6; DHHC-6; ZNF376; transmembrane protein H4; zinc finger DHHC domain-containing protein 6; zinc finger protein 376; zinc finger DHHC-type containing 6		

For research use only.

GENE DESCRIPTION	zinc finger, DHHC-type containing 6
GENE ID	64429
ORIGIN SPECIES	Human
ACCESSIONS	BC007213
VECTOR	pDONR223
NCBI SYMBOL	ZDHHC6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]