

## ABHD14B-abhydrolase domain containing 14B Gene

Reference: PTXBC007234

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | protein ABHD14B; CIB; HEL-S-299; CCG1-interacting factor B; abhydrolase domain-containing protein 14B; alpha/beta hydrolase domain-containing protein 14B; cell cycle gene 1-interacting factor B; epididymis secretory protein Li 299; abhydrolase domain cont |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | abhydrolase domain containing 14B                                                                 |
| <b>GENE ID</b>          | 84836                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC007234                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | ABHD14B                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggcagcaagcgtggagcagcgcgagggcaccatccagggtgcagggccagggccctcttctccgagagggccctgccggcagtgggcaggtcgtctctgtactgctgctgcatggtattcgttctcctccgagacctggcagaacctgggtacactgcacaggtggcccagg  
ctggctaccgggctgtggccattgacctgccaggtctggggcactccaaggaaagcagcagccctgcccctattggggagctggccctggcagcttctggcggtgtgtgtggtatgcttggagctgggcccccggtgtgatcagttccatcactgagtggtcatgactccctgccctt  
cttcacggccctggctccagctcccggtctgtgccaagtggcccatctgcactgacaaaatcaatgctgccaactatgccagtgtgaagactccagctctgattgtatgagaccaggacccatgggtcagaccagctttgagcacctgaagcagctgcccaccaccgggtg  
ctgatcatgaagggggcggggcaccctgttacctggacaaaccagagagtggtacatacagggtgctggaacttctgcaggggctccagtga