

C14orf128-chromosome 14 open reading frame 128 Gene

Reference: PTXBC007251

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf128; ARHGAP5 antisense RNA 1 (head to head)		
For research use only.			

GENE DESCRIPTION	chromosome 14 open reading frame 128
GENE ID	84837
ORIGIN SPECIES	Human
ACCESSIONS	BC007251
VECTOR	pDONR223
NCBI SYMBOL	C14orf128
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggccctgattcagtagctgtgaaagtggagcgagaagccaagacgtggatcgagaagccccgtggggcgggtctcaggggtggccaaaagacccctgttcacgcactaccagcctaacttaggaactgtcgtccactggctttattatactgccttaa