



Reference: PTXBC007276

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL-33; HEL-S-72p; HSC54; HSC71; HSP71; HSP73; HSPA10; LAP-1; LAP1; NIP71; heat shock cognate 71 kDa protein; LPS-associated protein 1; N-myristoyltransferase inhibitor protein 71; constitutive heat shock protein 70; epididymis luminal protein 33; epididy		

For research use only.

GENE DESCRIPTION	heat shock 70kDa protein 8
GENE ID	3312
ORIGIN SPECIES	Human
ACCESSIONS	BC007276
VECTOR	pDONR223
NCBI SYMBOL	HSPA8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr