

RPS11-ribosomal protein S11 Gene

Reference: PTXBC007283

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	S11; 40S ribosomal protein S11; ribosomal protein S11		
For research use only.			

GENE DESCRIPTION	ribosomal protein S11
GENE ID	6205
ORIGIN SPECIES	Human
ACCESSIONS	BC007283
VECTOR	pDONR223
NCBI SYMBOL	RPS11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggacattcagactgagcgtgcctacaaaagcagccgacctcttcaaaacaagaagggtcctgctgggagaaactggcaaggagaagctcccgctgactacaagaacatcggctcgggtcgaagacccaaggaggctattgaggcacctacattgacaagaat
gccccttcaactggtatgtgtccattcagggcgatcctctctggtggtgacaaagatgaagatgcagaggacctgtcatccgccgagactatctgactacatccgaagtacaaccgcttcgagaagcgccacaagaacatgtctgtacacctgtccccctgctcaggagcgtc
cagatcgggtgacatcgtcacagtggcgagtgccggcctctgagcaagacagtgcgttcaacgtgctcaagggtcaccaaggctgccggcaccaagaagcagttccagaagtctga