

MGC15634-hypothetical protein MGC15634 Gene

Reference: PTXBC007286

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein MGC15634
GENE ID	84841
ORIGIN SPECIES	Human
ACCESSIONS	BC007286
VECTOR	pDONR223
NCBI SYMBOL	MGC15634
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttgacctgcttccatgtctcaaacattagtcctacaagaaccctatgaggtttatctctttcacccattttgagagaagaaaagttcagctcagagagctgcaggaaattggggacatctccaagtcacagccaatagggtggcagtcagtggtgtcttgaagggaacaattg
agcttctgaacagctactccagaactatgggtgctgtggtgaagtcctggttaggagcctag