

PHF5A-PHD finger protein 5A Gene

Reference: PTXBC007321

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	INI; Rds3; SAP14b; SF3B7; SF3b14b; bK223H9.2; PHD finger-like domain-containing protein 5A; PHD finger-like domain protein 5A; PHD-finger 5a; splicing factor 3B associated 14 kDa protein; splicing factor 3b, subunit 7; PHD finger protein 5A		

For research use only.

GENE DESCRIPTION	PHD finger protein 5A
GENE ID	84844
ORIGIN SPECIES	Human
ACCESSIONS	BC007321
VECTOR	pDONR223
NCBI SYMBOL	PHF5A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctaacaatcatcctgattgatcttttgcgcaagcaggctgtgtgtccatcgaagactgtgtgaaaaatgtgatggcaagtgtgtgattgtgactcctatgtgcgtccctgcactctgggtgcgcataatgtgatgagtgtaactatggatcttaccagggcgctgtgtgatctgtggag
gacctgggtctctgatgcctattattgtaaggagtgaccatccaggagaaggacagagatggctgcccagaatgtcaatctggggagctctaagacagacctctctatgaacgcaaaaaatagcgttcaagaagaggta