

TMEM50A-transmembrane protein 50A Gene

Reference: PTXBC007341

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IFNRC; SMP1; transmembrane protein 50A; cervical cancer oncogene 9; small membrane protein 1		

For research use only.

GENE DESCRIPTION	transmembrane protein 50A
GENE ID	23585
ORIGIN SPECIES	Human
ACCESSIONS	BC007341
VECTOR	pDONR223
NCBI SYMBOL	TMEM50A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggatttctagaggcttgagatgctcagaatgcattgactgggggaaagcgcaatactattgctccattgctgctgggtgactattttacaggctgggtgattatagatgcagctgtattatcccatgaaagattcaaccactcattaccatgcctgtggtgttagcaacatagccttcctaatgattaatgcagatcgaatggacaagtcgaggtgatagttacagtgaagggtgctgggtcaaacagggtgctgcatttggcttttctgttttcattgttgcccttggatctctgattgcattatgtggattcttttggagggttatgttctaaagaaaaacatagtat accctggaattgctgtattttccagaatgccttcatttttggagggtggttttaagtttggccgactgaagacttatggcagtgga