

MGC16075-hypothetical protein MGC16075 Gene

Reference: PTXBC007354

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein MGC16075
GENE ID	84847
ORIGIN SPECIES	Human
ACCESSIONS	BC007354
VECTOR	pDONR223
NCBI SYMBOL	MGC16075
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagagaccctatgactcccgacatgcattgggtcaagtcacagtagaggaccactgagctgggcaccatgtggacctatagcaggcagcaggcagccctctcctttgaggttcctcctgcaactacgaccccgaaaaatggatccagtaa