



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

STX3-syntaxin 3 Gene

Reference: PTXBC007405

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

STX3A; syntaxin-3; syntaxin 3A; syntaxin 3

For research use only.

GENE DESCRIPTION	syntaxin 3
GENE ID	6809
ORIGIN SPECIES	Human
ACCESSIONS	BC007405
VECTOR	pDONR223
NCBI SYMBOL	STX3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaccgtctgagcagctgaagcccaagcagctgcacacaggaatgatatactgactcggttgattgctatcgacaacacgcctttatgacaggctctttctgagattgaggaaactcgcttaacattgacaagatctagaactgtagaggagctaaagaaactctacagt
atcattctctctgaccgattccagagccaaaacaaaggatgacgtagcagctcagcactgagtaagaagaaagggccaacacgtccggaacaaactgaagcattggaagacatattgaagaagatgaggtcaggtcatcgcgacacttcggattcgaaatcccgacactc
tgtctcttctcggaatttgttggagtgatgacacaaatacaatgaagctcaaatgactccgagacgcagcagaagggcggaalccagcgcgcagctcgaattactcgcaaaaagacaacgcgtgagggctgaggaagcattgttggagatgttggagatgttggcaaccggccatcttcaactctgggat
cattgacagacatttccaaagccctcagctgagattgaggtgagggagcagacagactgttgaggctgagagcagcatcgaagcgtcttcacatgctttatggacatcccatcggttggagaatcatggttgaatctgtagatacatagatgtgaatgtcatgcacacgtggac
cagctcgagaagccagagatgaaccgaaaaaacctctgaaatcacacagcctcagcgcccggaagaattgataattatcattctgctagtgtctgcttattgacattttagcatctatattgacttcttcggctcgaaatta