

NMB-neuromedin B Gene

Reference: PTXBC007407

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	neuromedin-B; Neuromedin-B-32; neuromedin-beta		
For research use only.			

GENE DESCRIPTION	neuromedin B
GENE ID	4828
ORIGIN SPECIES	Human
ACCESSIONS	BC007407
VECTOR	pDONR223
NCBI SYMBOL	NMB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccggcgggcgggcgctcggtggtggcagcctctgctcttcgacctgctgctgccggcgctgccccgctcagctgggatctccggagccccgcagccgagccagcaagatccgagtcgactcgcgaggcaacctctgggccaccggtcacttcattgggcaagaa
gagtcctggagcctccagccccatctggggacagctacccacacctccctgagggaccagcgactgcagctgagtcctgctgctggaatcctctgctaaagaaggctctggcgctgagcctcagcccccgcaccccaatccagaggctgctgtacaaatactgcag
aatgacaccaataatggggcagacacaacagcgtggcttagattgtgccccccaggggaaggctgctgaatgggacctgttgatggcccatctggatgtaaatctga