



Reference: PTXBC007413

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dJ1187M17.5; FAST kinase domain-containing protein 5, mitochondrial; FAST kinase domain-containing protein 5; FAST kinase domains 5		

For research use only.

GENE DESCRIPTION	FAST kinase domains 5
GENE ID	60493
ORIGIN SPECIES	Human
ACCESSIONS	BC007413
VECTOR	pDONR223
NCBI SYMBOL	FASTKD5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]