



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ACP1-acid phosphatase 1, soluble Gene

Reference: PTXBC007422

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

HAAP; LMW-PTP; LMWPTP; low molecular weight phosphotyrosine protein

phosphatase; LMW-PTPase; acid phosphatase of erythrocyte; adipocyte acid

phosphatase; cytoplasmic phosphotyrosyl protein phosphatase; low molecular weight

cytosolic acid phosphatase; Prot

For research use only.

GENE DESCRIPTION	acid phosphatase 1, soluble
GENE ID	52
ORIGIN SPECIES	Human
ACCESSIONS	BC007422
VECTOR	pDONR223
NCBI SYMBOL	ACP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggaacagcgctaccaagctccgtgctgttggctgtgtcggtaaacatttgcgtaaccctatgcagaagcagtttcaggaaactgtgaaccgatcaaaacatctcagagaattgggtcattgacagcggtgctgctgtttctgactggaaactggggccgggtcccgacaccaagagctgtga
gctgcctcaaaagaatctgactctacacagcccttaaaagaagaagactgtccacattgattatactatgattgagatgaagacactgcagagattgaatgaagaaataacagataaacctcgaagactctaaattgaactctgggagctatgactcaacaaa