

TMEM101-transmembrane protein 101 Gene

Reference: PTXBC007438

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 101		
For research use only.			

GENE DESCRIPTION	transmembrane protein 101
GENE ID	84336
ORIGIN SPECIES	Human
ACCESSIONS	BC007438
VECTOR	pDONR223
NCBI SYMBOL	TMEM101
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtcgaagataggttcgagacgggtgatgttcagctgatcatgagttgggttcggtgctgcacacgtgccccttttgggctgcttcagccagctcatgtgtacgtgagagggctgagggcagccgggaagcccgacatccagtgcttacctgtatttcacatggggc
agccgtgctgtgcgcatgtttcatgtcctttggcgtgaagcggcgtggttcgctggggccgcactccaattggcattagcacctacgcgcctacatcggggctacgtccactacgggactggctgaaggtcgtatgtactcgcgcacagttgccatcatcggcgtttctgt
gttggccagcgggtgctggggagctgtaccgcggaaacctcgcagccgctcctcagttccaccggccaggtgttctgtggtatctactcatctgtgtggcctactcactgcagcacagcaaggagaccggctggcgtatctgaacctctccaggaggaggctgatgatccagc
tgttcttctgtgtatggcactcggccctggcctttctgtcaggtactactgtgacctcgtgcccagatcctggctgtactgctgccccctgtcatgctgtcattgatggcaatgttgcctactggcacaacacgcggcgtgttgagttctggaaccagatgaagctccttgagagagtg
gggcactctcgaactgctgtcatcctggcactgatggctga