

WBP2-WW domain binding protein 2 Gene

Reference: PTXBC007452

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GRAMD6; WBP-2; WW domain-binding protein 2; WW domain binding protein 2		
For research use only.			

GENE DESCRIPTION	WW domain binding protein 2
GENE ID	23558
ORIGIN SPECIES	Human
ACCESSIONS	BC007452
VECTOR	pDONR223
NCBI SYMBOL	WBP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctcaacaagaatcactcggaggcgaggagtgatcgtcaataacaccgagagcctcctaatgtcctatgatcacgtggaactcacattcaatgacatgaagaacgtgccagaagcctcaaggaccaaaggaagcactgtctaccttacccttaccgggtcatcttctgtc
caagggaaggaatgcatgctcttcatgcatcttattatctcatgaagagctgtgagatcaagcagccgtatttgggtgcaaaactacatcaagggaacagtgaaaggcggaagcgagggtggtgggaaggctctgcttctacaagttgactttcacggcagggggcgccatgag
ttcgacagcggatgctccagggtggcatctcaagcctccagaggtgaagtcctcagtgagcctatggctactcttaccatgccagcggggcctatgtctatccccgccagtcgccaatggaatgtacccctgcccctggctacccctatccaccgccccacctgagttctatccagg
acccccatgatggacggggccatgggatactgagccccaccacggcctacccctggggccatgggaacctccggtcagcggcccgatgtccctccactctgcagccgaagccaagggccgcagaagcagccgcagcgcctattacaaccaggcaatctcacaacgtct
acatgccacgagccagccggccacctccctactaccacgggaagataagaagaccagtag