



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ZNF503-zinc finger protein 503 Gene

Reference: PTXBC007494

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

NOLZ-1; NOLZ1; Nlz2; zinc finger protein 503

For research use only.

GENE DESCRIPTION	zinc finger protein 503
GENE ID	84858
ORIGIN SPECIES	Human
ACCESSIONS	BC007494
VECTOR	pDONR223
NCBI SYMBOL	ZNF503
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcacagcgccctcgcttctgccctaagaagcagtaagacacagcggcgccggcgccggcgaggcgggagcgccggcggtgcagacctgcctggacacagcgcgctctctggaataatgcctccggccccggccaggctcgtctcccgccggcgagcaccagacctttgtgcac
gccgtgccccctctgacctctgcccagcgccaacgctctgccaatcaaggtgctgaagatgctgacggcgacgaactggcccaattttgcacccgtagtactgcagccccctgccttccacgcccgtgacccccatcagctcgtatgccaaagagagcccgctggcgctgttgcgcg
aaacatgttcgcagatcgggaagcccgacacctgcctccctgcgcggccacacccctctaccctcagcgctttatgctctccaaacgacctccccacatctgcaactgggtgtcggcccaagggcgctgcgcgaacggctccacgtctccacgtctcgaagagctgcgtgacccactgcg
accatacgcgaattccggcagacaaactgctctggctcaccacagctctgctctgcgcagcggcgccggcgacagctgttcacatgacacatccccactggcgcaacggggagccctggagcgtgcgcgtgcgcagccccaccacgctgcgtgagactcag
cagcgtctaccacccctctcaagaagcgcctctccagcgtctgcgcctctccagctgcgcggcgccgacccggcgactactactccctcctacgctccctacgacagacagagatgaccacgcctcgccgctcgtggtatcagta