

RPS20-ribosomal protein S20 Gene

Reference: PTXBC007507

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	S20; 40S ribosomal protein S20; ribosomal protein S20		
For research use only.			

GENE DESCRIPTION	ribosomal protein S20
GENE ID	6224
ORIGIN SPECIES	Human
ACCESSIONS	BC007507
VECTOR	pDONR223
NCBI SYMBOL	RPS20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttttaaggataccggaaaaacacccgtggagccggaggtggcaattcacgaattcgaatcacccctaacaagccgcaacgtaaaatccttgaaaaaggtgtgtgctgactgataagaggcgcaaaagaaaagaatctcaaaagtgaaggaccagttcgaatgcctaccaagact
ttgagaatcactacaagaaaaactccttgggtgaagggttctaagacgtgggatcggttccagatgagaattcacaagcgactcattgacttgcacagtcctctgagattgttaagcagattacttccatcagttgagccaggagttgaggtggaagtcaccattgcagatgcttaa