

SPAG9-sperm associated antigen 9 Gene

Reference: PTXBC007524

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT89; HLC-6; HLC4; HLC6; JIP-4; JIP4; JLP; PHET; PIG6; C-Jun-amino-terminal kinase-interacting protein 4; JNK interacting protein; JNK/SAPK-associated protein; Max-binding protein; c-Jun NH2-terminal kinase-associated leucine zipper protein; cancer/testis		

For research use only.

GENE DESCRIPTION	sperm associated antigen 9
GENE ID	9043
ORIGIN SPECIES	Human
ACCESSIONS	BC007524
VECTOR	pDONR223
NCBI SYMBOL	SPAG9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtcctggatgcctgctgtttgtgtttgctttgttggcggggcggtggtcattatctgctatcttagtactctctgtttgtctgttgacattttctatttctaccggtgtgccagctctgacgcagaacctaccaaggatactcagaaaaaacgccctatatcattaggaatttcc
cattacctgctggagatggattgcttacctgatctcagaaaggaggagagaccctggatcagcaatggaaattcaggaattaaagtaaccacgttctcaccagcctgaaggtcagcaatagctcctgaacctcagaaggctgtagaacaggaggtgagaatggtgtcttctaaca
ttttgcaaaaagtatactaa