



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KIAA1984-KIAA1984 Gene

Reference: PTXBC007542

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

KIAA1984; PARF; bA216L13.7; coiled-coil domain-containing protein 183; coiled-coil domain containing 183

For research use only.

GENE DESCRIPTION	KIAA1984
GENE ID	84960
ORIGIN SPECIES	Human
ACCESSIONS	BC007542
VECTOR	pDONR223
NCBI SYMBOL	KIAA1984
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]