



Reference: PTXBC007576

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DCAL; IRBIT; PPP1R78; PRO0233; XPVKONA; adenosylhomocysteinase 2; DC-expressed AHCY-like molecule; IP(3)Rs binding protein released with IP(3); S-adenosyl homocysteine hydrolase homolog; S-adenosyl-L-homocysteine hydrolase 2; S-adenosylhomocysteine hydrol		

For research use only.

GENE DESCRIPTION	S-adenosylhomocysteine hydrolase-like 1
GENE ID	10768
ORIGIN SPECIES	Human
ACCESSIONS	BC007576
VECTOR	pDONR223
NCBI SYMBOL	AHCYL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]