



Reference: PTXBC007596

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATG18B; Atg21; CGI-50; WIPI-2; WD repeat domain phosphoinositide-interacting protein 2; WD40 repeat protein interacting with phosphoinositides 2; WIPI49-like protein 2; WD repeat domain, phosphoinositide interacting 2		

For research use only.

GENE DESCRIPTION	WD repeat domain, phosphoinositide interacting 2
GENE ID	26100
ORIGIN SPECIES	Human
ACCESSIONS	BC007596
VECTOR	pDONR223
NCBI SYMBOL	WIP12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]