

PTMS-parathymosin Gene

Reference: PTXBC007616

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ParaT		
For research use only.			

GENE DESCRIPTION	parathymosin
GENE ID	5763
ORIGIN SPECIES	Human
ACCESSIONS	BC007616
VECTOR	pDONR223
NCBI SYMBOL	PTMS
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccaggagagagagtcctcaggaaaaggccccagggacctaaggagaagaaggagaagggtggaggagaaggcaagccggaaagagcgaagaagaagaagtgggtggagcctccagtccttgggtgatgggtggcaaggcatgcagccagcctgaggctactccctc
tcctgtgtccccacaggaggaggagaacggggctgaggaggagaagaagaactgccgaggatggagaggaggagaagatgaagggaagaagaagatgaggagaagaagaaggatgatgaaggggccgcgtga