



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

DKK3-dickkopf homolog 3 (*Xenopus laevis*) Gene

Reference: PTXBC007660

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	REIC; RIG; dickkopf-related protein 3; RIG-like 5-6; RIG-like 7-1; dickkopf 3 homolog; dickkopf homolog 3; dickkopf-3; dkk-3; hDkk-3; regulated in glioma; dickkopf WNT signaling pathway inhibitor 3		

For research use only.

GENE DESCRIPTION	dickkopf homolog 3 (<i>Xenopus laevis</i>)
GENE ID	27122
ORIGIN SPECIES	Human
ACCESSIONS	BC007660
VECTOR	pDONR223
NCBI SYMBOL	DKK3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr