

NDUFB9-NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 9, 22kDa Gene

Reference: PTXBC007672

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B22; CI-B22; LYRM3; UQOR22; NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 9; LYR motif-containing protein 3; NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 9, 22kDa; NADH-ubiquinone oxidoreductase B22 subunit; complex I B22 subunit; NADH:u		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 9, 22kDa
GENE ID	4715
ORIGIN SPECIES	Human
ACCESSIONS	BC007672
VECTOR	pDONR223
NCBI SYMBOL	NDUFB9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgttcttggcgtcggaccctacctgacccatcagcaaaagggtgttgcggctttataagcggcgctacgccacctcgagtcgttgcgtccagagagacaaataccgatacttgccttggatgagagcccggttgaagaacataagaatgaaaggatatggcgaaggccac
ccagctgctgaaggaggccgaggaagaattctgtaccgtcagcatccacagccatacatctccctgactctcctggggccacctcctatgagagatacgaattgctacaagggtccagaatgggtcttagatgactggcatcctcttgagaaggcaatgtatcctgattactttccaagaga
gaacagtggagaactgcgaggggaagctgggaacgagaggttaagcagctgcaggaggaaacccacctgtgtgctcttaactgaagcttgcctccctgcccgaagggaaggtgatttgcctccactgtgtgtgatttggaccagacccgggagcgcccatgtag