

ZNF324-zinc finger protein 324 Gene

Reference: PTXBC007717

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZNF324A; zinc finger protein 324A; zinc finger protein ZF5128; zinc finger protein 324		

For research use only.

GENE DESCRIPTION	zinc finger protein 324
GENE ID	25799
ORIGIN SPECIES	Human
ACCESSIONS	BC007717
VECTOR	pDONR223
NCBI SYMBOL	ZNF324
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctttgaggatgtggctgtgtacttctccaggaggagtggtgggtctctggacacagccagagggtccctgtaccgctgtgatgtagacaacttgcacttgtggctcgtgggactctccacctctcgacctcgtgtggtcatccaactggagcgtggcaggagccctgg
gttcccaagtggaaaggacacaacctgtccaggaccacctacaggaggcgcaacctgtgtctgtgagttgacagaggatagagatgtttctggagaatggccacgagcttccagataccccacctgggatgactactagcgtcttccctgttgcgggtgcctgccacagtgtaaaaa
gcctgcagagacaacgggtgtcctcccatctcgggagagaaaaccacgggggtgtcgtgtgaltacttgggagagctcctgttaggtctcaggcagtggtggcaagccagcgtcagcctgcgactgacctccccgttaggctccccaggggcgtccggcttagggaaaagacactc
acagagcatgcgttctggggaggcagccagggagcctgagcggcagaacctgtgcacagagggtccctgggagaaccttgggagcgccagcagctggaggctgccggcggtcggggacatcaccgaatgggtgcagtttggcaggagcctcatagactcctcgtggcc
aggagccctcgactgggacgagctggcgaggctcttcacgtggggagaaagtccttcgaatgcaggcgctgcagaaaagtgctgtgaagagctccgacctctcaagcactacgcacccgtgtgcggcaaggccttccggcatag