

COX6A1-cytochrome c oxidase subunit VIa polypeptide 1 Gene

Reference: PTXBC007723

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMTRID; COX6A; COX6AL; cytochrome c oxidase subunit 6A1, mitochondrial; COX VIa-L; cytochrome C oxidase subunit VIa homolog; cytochrome c oxidase polypeptide VIa-liver; cytochrome c oxidase subunit VIA-liver; cytochrome c oxidase subunit VIa polypeptide 1		

For research use only.

GENE DESCRIPTION	cytochrome c oxidase subunit VIa polypeptide 1
GENE ID	1337
ORIGIN SPECIES	Human
ACCESSIONS	BC007723
VECTOR	pDONR223
NCBI SYMBOL	COX6A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggtagttgggtgtctcgtttctcggctgctgggtcgggtccgccacagctggggcgccctatgtcagtggtggcccatggcgaagagggctcagctcgcgtggaagactctcacccttctcgcgcgtcccccgggtggcagtcagcatgctgaatgtgtacctgaagt
cgccaccggagagacgagagacccgagttcatgcctaccgccatctccgcacagaccagccgtttccctggggagatggttaaccatactctattccataaccctcatgtgaatccacttccaactggctacgaagatgaataa