



Reference: PTXBC007741

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DLK; DLK-1; Delta1; FA1; PREF1; Pref-1; ZOG; pG2; protein delta homolog 1; delta-like 1 homolog; fetal antigen 1; preadipocyte factor 1; securedeltin; delta like non-canonical Notch ligand 1		

For research use only.

GENE DESCRIPTION	delta-like 1 homolog (Drosophila)
GENE ID	8788
ORIGIN SPECIES	Human
ACCESSIONS	BC007741
VECTOR	pDONR223
NCBI SYMBOL	DLK1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]