

CKS1B-CDC28 protein kinase regulatory subunit 1B Gene

Reference: PTXBC007751

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CKS1; PNAS-16; PNAS-18; ckshs1; cyclin-dependent kinases regulatory subunit 1; CDC2-associated protein CKS1; CDC28 protein kinase 1; CDC28 protein kinase 1B; CKS-1; NB4 apoptosis/differentiation related protein; PNAS-143; cell division control protein CKS		

For research use only.

GENE DESCRIPTION	CDC28 protein kinase regulatory subunit 1B
GENE ID	1163
ORIGIN SPECIES	Human
ACCESSIONS	BC007751
VECTOR	pDONR223
NCBI SYMBOL	CKS1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgacacaacaatttactattcggacaaatagacgacgaggagtttgatcgcacatgtcatgtctgccaaggacatagccaagctgtgtcctaaacccatctgatgtctgaatctgaatggaggaatcttggcgttcagcagatcaggatgggtccattatgatccatgaac
cagaacctcacatcttctgttccggcgccactaccgaagaaccaagaagaatga