



Reference: PTXBC007800

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DFNB94; SLM5; asnRS; asparagine tRNA ligase 2, mitochondrial (putative); deafness, autosomal recessive 94; asparaginyl-tRNA synthetase 2, mitochondrial (putative)		

For research use only.

GENE DESCRIPTION	asparaginyl-tRNA synthetase 2, mitochondrial (putative)
GENE ID	79731
ORIGIN SPECIES	Human
ACCESSIONS	BC007800
VECTOR	pDONR223
NCBI SYMBOL	NARS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1