

GPN2-GPN-loop GTPase 2 Gene

Reference: PTXBC007815

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATPBD1B; GPN-loop GTPase 2; ATP-binding domain 1 family member B		
For research use only.			

GENE DESCRIPTION	GPN-loop GTPase 2
GENE ID	54707
ORIGIN SPECIES	Human
ACCESSIONS	BC007815
VECTOR	pDONR223
NCBI SYMBOL	GPN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaggggcccgcaccacggccttcgggcaggcggatgacggcccgccgggctcagggaagaccacgtactgcctgggcatgagtgagttcctgcgcgcgtggccggcgcggtgggtgaacctggaccggccaacgaggggctgccgtacgagtgccgtg
gacgtggcgagctggtgggctggcgacgtgatggacgctgcgctggggcccaacggcgccctgctctactgcatggagtaacctggagccaacctggactggctgctgccaagctcgacccctccggccactactctcttgcactgcccaggcagtgagctc
tgcacgcatcagggcgcttgcgcagcatcttcccaaatggcgagtgaggacacagggctgactgcccgtccacctcgtggattctcactactgcacagacctgccaagtcaattcagtagtctgttacctccctggccaccatgctgcacgtggaaactgccccacatcaacctcttcca
agatggacctcattgagcattatgggaagctggccttcaacctggactactacacagaggttctgacgtctcctacgtcctgaccacctgggtctgacctttcttccgactaccgccagctcaatgagaagctagtcagctcctcgaagactatagcctgtctctttatcccttcaa
catccaggacaaggagagcatccagcgagtcctgagcgtgtggataaagccaatggatactgtttcggagcccaagagcagcgaagcttggagccatgatgtctgcccgaatggagccgacttccattctctccacactgggcatccaggagaagctacgtggcaccctcgaacc
agtcagtgagcaggaagccatgcagctgtag