

C3orf34-chromosome 3 open reading frame 34 Gene

Reference: PTXBC007827

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf34; MOSPGF; centrosomal protein of 19 kDa; centrosomal protein 19kDa; centrosomal protein 19		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 34
GENE ID	84984
ORIGIN SPECIES	Human
ACCESSIONS	BC007827
VECTOR	pDONR223
NCBI SYMBOL	C3orf34
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgtgcactgccaagaaatgtggattagttcagcctccagctattatcttaatctatgagagtgaatcaagggaattcgccagcgattatgccagttcgaactttcaaagtttcagattgcaccagagctgctgaacaattaaagaataatccgcgacacaagagttacctag
aacaagtatccctgagcgagctagagaagctatcagttttttacgaggttactgtcggggcagagtcgagaaacaatggaacaaattcaacgggaacaaccattgatcctgaggaagacctgaacaaactagatgacaaggagcttgccaaaagaaaagagcatcatggatgaact
ttttgagaaaaatcagaagaagaaggatgatccaaattttgtttatgacattgaggttgaatttcacaggacgatcaactgcagtcctgtggctgggacacagagtcagctgatgagttctga