

YIPF5-Yip1 domain family, member 5 Gene

Reference: PTXBC007829

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein YIPF5; FinGER5; SB140; SMAP-5; SMAP5; YIP1A; YIP1 family member 5; YPT-interacting protein 1 A; five-pass transmembrane protein localizing in the Golgi apparatus and the endoplasmic reticulum 5; golgi membrane protein SB140; smooth muscle cell ass		

For research use only.

GENE DESCRIPTION	Yip1 domain family, member 5
GENE ID	81555
ORIGIN SPECIES	Human
ACCESSIONS	BC007829
VECTOR	pDONR223
NCBI SYMBOL	YIPF5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcaggcttgaacacggaattctaccagacaagttacagcatcagtcacagcagtcctatgattatggaggagtgaggaccctatagcaaacagtagctggctatgactatcgcagcaaggcagattgtccctccagacatgatgcagccacaacagccat
acaccgggcagatttaccagccaactcagcgcatactcagcttcacctcagcctttctatggaaacaactttgaggatgagccacctttattagaagagttaggtatcaattttgaccacatctggcaaaaaacactaacagtattacatccgttaaaagtagcagatggcagcatcatgaatg
aaactgatttggcaggtccaatgggttttgccttgcctttggagccacattgctactggctggcaaaatccagtttggctatgtatcgggacagtcgaattggatgctaggaatgtttttatataacttaagtagtatgacaggtgtttcatttgggtgtggtggaagtgcttggatattgtcttc
tgccatgacatcttccagcttgcagtgatatttttgcagggaatgttaggaatcttctactgctgggattatggatggtgtatgtttctgctccaaaatattttctgcatagccatggaaggacagcaacttttagcatatccttgcgcttggatattgagcttgcctgat
ttccgtctttga