

PDE6D-phosphodiesterase 6D, cGMP-specific, rod, delta Gene

Reference: PTXBC007831

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	JBTS22; PDED; retinal rod rhodopsin-sensitive cGMP 3',5'-cyclic phosphodiesterase subunit delta; GMP-PDE delta; phosphodiesterase 6D, cGMP-specific, rod, delta; protein p17; phosphodiesterase 6D		

For research use only.

GENE DESCRIPTION	phosphodiesterase 6D, cGMP-specific, rod, delta
GENE ID	5147
ORIGIN SPECIES	Human
ACCESSIONS	BC007831
VECTOR	pDONR223
NCBI SYMBOL	PDE6D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcagccaaggacgagcgggcccaggagatcctgagggcctcaactaaattggatgaaccttcgggatgctgagacagggaagatactctggcaagggaacagaagacctgtctgtccctggtgtggagcatgaagcccggtgtcccaagaaatcctcaagtgcaggcagtgctc
tcgagaacttaattttctcgacagaacaatggaaaaatccgctgggaacaaaagtcttcaaaaggcgaatgcctagaagaatggttcttcgagtttgctttgtgatccctaactccacaataacctggcagtcctgatagaggcagcaccgagtcacagatgatgccagcaagcg
tcttaactgggaacgttatcatagaacaagaatttttgacgacgatctcttgaagcacatccagatgagactttctatgtttga