

FBXO44-F-box protein 44 Gene

Reference: PTXBC007832

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FBG3; FBX30; FBX6A; Fbx44; Fbxo6a; F-box only protein 44; F-box gene 3; F-box protein FBX30; F-box/G-domain protein 3; F-box protein 44		

For research use only.

GENE DESCRIPTION	F-box protein 44
GENE ID	93611
ORIGIN SPECIES	Human
ACCESSIONS	BC007832
VECTOR	pDONR223
NCBI SYMBOL	FBXO44
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtggggaacatcaacgagctgcccagaaacatctgctggagctgttcacgacgtgcccggccgagctgctgaactgccgctgtctgcagcctctggcgggacatcagacctctgacccctctggaacgcaagtgcctgcgagagggctcatcactgagga
ctgggaccagcccggtggccgactggaagatcttacttcttacttcttactgagcctgcacaggaacccctgcacaaacccgtgcgctgaagaggggttcgagttctggagcctggatgtgaatggagcgaatgagtggaaggtggaggtatctctcagaccagaggaaggaattcccaatg
accaggttcgagccagccagattgcgggtccaagtaccagctgtgcttcagctcctgtcgtccgcgcacgcgcctctgggacatccagccagaccggcgacatccagcagaagagcgaatgccaaagtgaggagggtctccacattctccaactaccggccggcgtc
cgctacatctgggttcagcagcgccgctggacactcattactgggcccgtgttacggccggaggtcaccacagcagcaccatccggcccccgtgccctgacacccctgagccccatctgctgaacctga